



Surfaces & Pads

Allen® Perioperative Surfaces & Pads

Allen Medical offers a variety of pads to fit your standard OR tables, armboards and stretchers. Allen has over 10 years of experience in the area of pressure management. The engineered foams in our stretcher and OR table pads were cutting edge when first developed and our products have remained at the forefront of technological innovation ever since - with new radiolucent foam formulations and multi-density pad layouts.


Why OR Table Surfaces Matter...

Healthcare professionals are implementing new prevention initiatives and standards for pressure ulcer prevention.


From a surge in recent research pertaining to pressure ulcers to the Medicare Guidelines that took effect in October 2008, the issue of hospital-acquired pressure ulcers is a popular topic for many healthcare professionals. Additionally, professional organizations like AORN and NPUAP are critical to advancing awareness of the issue through their release of new standards and practices to help prevent and treat hospital-acquired pressure ulcers.

- Static surfaces increase surface support areas and decrease pressures and shear forces.⁴
- Old-fashioned "egg crate" mattresses offer no advantage.⁴
- A study of 1,128 surgical patients, revealed the incidence of nosocomial pressure ulcers related to the operating room (OR) was 23%.⁵
- A multi-layer pressure redistribution surface was found to be statistically significant in reducing pressure ulcer incidence when compared with the "Standard" OR table pad.⁵
- Surgical patients account for 34-42% of all hospital-acquired pressure ulcers.²
- In an experiment involving 323 patients between 2000 and 2004, 38% of patients who had the standard pad during surgery developed the sores. But only 7% of patients who had the pressure-relieving pad developed pressure ulcers.⁷
- Various studies have suggested that pressure ulcers may occur within a 1 to 4 hour time frame, well within the range of current surgical procedures.⁶

Pressure Ulcer Prevalence



Cost Of OR Acquired Pressure Ulcers²



Hospital Acquired Pressure Ulcers²

1. Aronovitch SA. Intraoperatively acquired pressure ulcer prevalence: a national study. *J Wound Ostomy Continence Nurs.* 1999;26(3):130-136.

2. Beckrich K, Aronovitch SA. Hospital-acquired pressure ulcers: a comparison of costs in medical vs. surgical patients. *Nurs Econ.* 1999;17(5):263-71.

3. Stotts N. Predicting pressure ulcer development in surgical patients. *Heart and Lung.* 1988;17(6):641-6.

4. Pressure Ulcers (Pressure sores; Bedsores; Decubitus Ulcers; Decubiti). *The Merck Manuals, Online Medical Library.* Available at: <http://www.merck.com/mmpe/sec10/ch126/ch126a.html> (accessed Sept. 2015).


5. Lummus AC. Perioperative Pressure Ulcer Assessment and Prevention: Efficacy Study of a Multi-Layer Pressure Relief Pad in the Operating Room (OR). University of Tennessee Health Science Center. 2004.

6. Baron S, MacFarlane G. Reducing Pressure Ulcer Risk In The Operating Room. White Paper. 2009. Allen Medical Systems, D-770444-A2.

7. Scott-Williams S. Perioperative Pressure Ulcer Prevention Program (PPUPP): An innovative effort to prevent pressure ulcers in surgical patients. *American Academy of Nursing.* <http://www.aannet.org/edge-runners-perioperative-pressure-ulcer-prevention-program> (accessed Sept. 2015).

Allen® Surfaces Technology

No Single Foam Provides Effective Pressure Management...



Allen Uses 4 Foams For Effective Pressure Management

Optimal pressure management requires a pad that is effective in four main support functions: 1-conforming to the patient's body, 2-providing proper support, 3-prohibiting bony prominences from bottoming out and 4-evening out irregularities and gaps in the table surface.

No single foam can support all four functions effectively. Listed below are the 4 foams Allen has proven effective to provide optimal pressure management in a single pad.

"Neither the position of the patient, nor the time spent on the OR bed can be altered easily during surgery."

Varying the type of pressure-reducing mattress may be the most appropriate method to prevent pressure ulcers during the perioperative period."

- Alexander's Care of the Patient in Surgery, 13th Edition



1
Memory Foam
Conforming layer



2
Low Density Elastic Foam
Support layer



3
High Density Elastic Foam
Anti-bottoming layer



4
Closed Cell Foam
Diaphragm layer


Allen Perioperative Surfaces: Pressure Maps

Model: Female
Weight: 117 lbs (53 kg)
Height: 63" (160 cm)


Alto




Astrus Elite



Deluxe




Standard



Model: Male
Weight: 180 lbs (82 kg)
Height: 68" (173 cm)


Alto



Astrus Elite



Deluxe




Standard




Model: Male
Weight: 300 lbs (136 kg)
Height: 73" (185 cm)


Alto




Astrus Elite



Deluxe



Standard



The Allen® Pressure Management Solution

Allen Pads Contain Multiple Foam Layers For Optimal Pressure Management



Alto Cross-Section

Alto™ Surface

Excellent Surface for Most Surgical Positions The Alto Surface is designed to offer superior pressure management in most surgical positions. Engineered, multi-layered and fully contoured, this surface allows bony prominences or “landing gear” to be protected in the Pressure Reducing Zones while other parts of the body are supported and cradled. This support system decreases the high pressure that bony prominences usually experience. The foot section is constructed with calf support and a special soft zone for heels. The anti-shear cover helps protect the patient’s skin.




Astrus Elite Cross-Section

Astrus Elite™ Surface

Effective Pressure Management Our Astrus Elite Surface has a 4-layer foam composite construction that protects patients in four ways. A top layer of memory foam conforms to the patient while the underlying foam layers support the patient’s weight and prevent bottoming out. Both the Astrus Elite and Astrus Elite Silver surfaces contain the same 4-layer advanced pad technology.

Allen Surfaces Selector Chart


		Foam Layers					Benefits				Pad Thickness							
		Specialty Layer	Conforming Layer	Support Layer	Anti-bottoming Layer	Diaphragm Layer	Velcro® Kit Included	Improved Pressure Management	Recommended for Long Surgeries	Anti-Shear Cover	2" (5 cm) Thickness	3" (7.6 cm) Thickness	3.25" (8.2 cm) Thickness	3.5" (9 cm) Thickness	4" (10 cm) Thickness	4.25" (10.8 cm) Thickness	5" (10 cm) Thickness	5.25" (13.3 cm) Thickness
Alto	OR Table Surfaces	●	●	●		●	●	●	●						●			
	Armboard Surfaces	●	●	●		●		●	●						●			
Astrus Elite & Silver	OR Table Surfaces		●	●	●	●	●	●	●		●			●		●		
	Armboard Surfaces		●	●	●	●		●	●		●	●		●	●	●	●	●
	Stretcher Surfaces		●	●	●	●		●	●				●					
Standard & Deluxe	OR Table Surfaces			●			●											
	Armboard Surfaces			●						●	●			●				
	Stretcher Surfaces			●						●			●					



Alto™ Surface

- Contoured to provide excellent pressure relief for all surgical positions
- Offered with integral gel layer (gray bottom) or without (black bottom)
- Velcro® kit included
- Non-conductive
- 4" Thickness

Alto OR Table Surface: Pressure Mapping Results¹




Clinical Data For Pressure Reducing Table Surfaces²


The Alto Surface outperformed standard OR table pads in a clinical study of 323 patients conducted by Suzy Scott-Williams, RN (study available upon request). The study group using the Alto Surface presented with dramatically fewer ulcers (7% ulcer rate) compared to the control group (38% ulcer rate).

Suzy Scott-Williams, RN, MSN, CWOCN, VA Medical Center, Memphis, Allan C. Lummus, PhD, Department of Preventive Medicine University of TN Health Science Center, Memphis. 2004

Control Group
(Standard Surface, n=176)



Study Group
(Alto Surface, n=147)




1. 6. Baron S, MacFarlane G. Reducing Pressure Ulcer Risk In The Operating Room. White Paper. 2009. Allen Medical Systems, D-770444-A2.
2. Study conducted by Suzy Scott Williams. Data concerning Alto Surface extracted by Allen Medical Systems, Inc. for development of this clinical sheet. The extracted data was reviewed and approved by Suzy Scott Williams and Allan Lummus for distribution by Allen Medical Systems, Inc.



Alto™ OR Table Pad Sets


	Alto Order Info 4" Pad	Alto Gel Order Info 4.25" Pad
OR Table Manufacturer / Model Name	Product Number	Product Number
Berchtold® B810, B710	H-006-12-100	H-016-12-100
Berchtold Operon® D850, D750 (long foot)	H-006-12-120	H-016-12-120
Berchtold Operon D850, D750 (short foot)	H-006-12-130	H-016-12-130
Maquet® 1130	H-006-10-110	H-016-10-110
Maquet 1130 Split Leg	H-006-10-100	H-016-10-100
Maquet 1130 Orthopedic	H-006-10-130	H-016-10-130
Maquet Alphamaxx® (long torso)	H-006-10-150	H-016-10-150
Maquet Alphamaxx (short torso)	H-006-10-160	H-016-10-160
Maquet Alphastar 1132	H-006-10-120	H-016-10-120
Midmark® 7100	H-006-06-100	H-016-06-100
Midmark Chic CST	H-006-06-120	H-016-06-120
Midmark Montreal	H-006-06-130	H-016-06-130
Getinge™ / Champagne® 4900B, 4900E	H-006-04-100	H-016-04-100
Getinge / Champagne 5100B, 5100E	H-006-04-130	H-016-04-130
Siemens® Universal	H-006-13-100	H-016-13-100
Skytron® 5001, 6001, 6500, 6600	H-006-05-130	H-016-05-130
Skytron 3100	H-006-05-100	H-016-05-100
Skytron 3500	H-006-05-120	H-016-05-120
Skytron 3600	H-006-05-190	H-016-05-190
Skytron 6700, 6701	H-006-05-200	H-016-05-200
Steris-Amsco 1080, 2080	H-006-03-100	H-016-03-100
Steris-Amsco 3080, 3085	H-006-03-110	H-016-03-110
Steris® Cmax® / Steris 4085	H-006-03-150	H-016-03-150
Steris SurgiGraphic® 6000	H-006-03-170	H-016-03-170

Individual pads for these sets are available. Contact Sales Dept. for pad part numbers.



Astrus Elite™ Surfaces

Effective Pressure Management Our Astrus Elite and Astrus Elite Silver™ Surfaces have a 4-layer foam composite construction that helps protect patients four ways. A top layer of conforming foam molds to your patients' anatomy. The pressure displacement foam middle layer transfers pressure off the bony prominences and the closed cell foam bottom layer prevents bottoming out. The anti-shear cover helps protect patient's skin. Offered with integral gel layer (gray bottom) or without (black bottom.) Available in 3" (7.6 cm) and 4" (10 cm). Velcro® kit included. Non-conductive.




Astrus Elite Silver Surfaces


Astrus Elite Silver Surface offers the same superior pressure management as the Astrus Elite Surface but now with silver embedded in the cover of the pad.

See page 128 for order information


Astrus Elite & Astrus Elite Silver Surface: Pressure Mapping Results




Astrus Elite
Supine Position



Standard Surface
Supine Position




Astrus Elite
Lateral Position



Standard Surface
Lateral Position

Model: Female
Weight: 205 lbs
Height: 66"





Astrus Elite™ OR Table Pad Sets

	Astrus Elite - Order Info			Astrus Elite Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
OR Table Manufacturer / Model Name	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Berchtold B810, B710	H-303-12-100	H-403-12-100	H-503-12-100	H-313-12-100	H-413-12-100	H-513-12-100
Berchtold Operon D850, D750 (long foot)	H-303-12-120	H-403-12-120	H-503-12-120	H-313-12-120	H-413-12-120	H-513-12-120
Berchtold Operon D850, D750 (short foot)	H-303-12-130	H-403-12-130	H-503-12-130	H-313-12-130	H-413-12-130	H-513-12-130
Maquet 1130	H-303-10-110	H-403-10-110	H-503-10-110	H-313-10-110	H-413-10-110	H-513-10-110
Maquet 1130 Split Leg	H-303-10-100	H-403-10-100	H-503-10-100	H-313-10-100	H-413-10-100	H-513-10-100
Maquet 1130 Orthopedic	H-303-10-130	H-403-10-130	H-503-10-130	H-313-10-130	H-413-10-130	H-513-10-130
Maquet Alphamaxx (long torso)	H-303-10-150	H-403-10-150	H-503-10-150	H-313-10-150	H-413-10-150	H-513-10-150
Maquet Alphamaxx (short torso)	H-303-10-160	H-403-10-160	H-503-10-160	H-313-10-160	H-413-10-160	H-513-10-160
Maquet Alphastar 1132	H-303-10-120	H-403-10-120	H-503-10-120	H-313-10-120	H-413-10-120	H-513-10-120
Midmark 7100	H-303-06-100	H-403-06-100	H-503-06-100	H-313-06-100	H-413-06-100	H-513-06-100
Midmark Chic CST	H-303-06-120	H-403-06-120	H-503-06-120	H-313-06-120	H-413-06-120	H-513-06-120
Midmark Montreal	H-303-06-130	H-403-06-130	H-503-06-130	H-313-06-130	H-413-06-130	H-513-06-130
Getinge / Champagne 4900B, 4900E	H-303-04-100	H-403-04-100	H-503-04-100	H-313-04-100	H-413-04-100	H-513-04-100
Getinge / Champagne 5100B, 5100E	H-303-04-130	H-403-04-130	H-503-04-130	H-313-04-130	H-413-04-130	H-513-04-130
Siemens Universal	H-303-13-100	H-403-13-100	H-503-13-100	H-313-13-100	H-413-13-100	H-513-13-100
Skytron 5001, 6001, 6500, 6600	H-303-05-130	H-403-05-130	H-503-05-130	H-313-05-130	H-413-05-130	H-513-05-130
Skytron 3100	H-303-05-100	H-403-05-100	H-503-05-100	H-313-05-100	H-413-05-100	H-513-05-100
Skytron 3500	H-303-05-120	H-403-05-120	H-503-05-120	H-313-05-120	H-413-05-120	H-513-05-120
Skytron 3600	H-303-05-190	H-403-05-190	H-503-05-190	H-313-05-190	H-413-05-190	H-513-05-190
Skytron 6700, 6701	H-303-05-200	H-403-05-200	H-503-05-200	H-313-05-200	H-413-05-200	H-513-05-200
Steris-Amsco 1080, 2080	H-303-03-100	H-403-03-100	H-503-03-100	H-313-03-100	H-413-03-100	H-513-03-100
Steris-Amsco 3080, 3085	H-303-03-110	H-403-03-110	H-503-03-110	H-313-03-110	H-413-03-110	H-513-03-110
Steris Cmax / Steris 4085	H-303-03-150	H-403-03-150	H-503-03-150	H-313-03-150	H-413-03-150	H-513-03-150
Steris SurgiGraphic 6000	H-303-03-170	H-403-03-170	H-503-03-170	H-313-03-170	H-413-03-170	H-513-03-170

Individual pads for these sets are available. Contact Sales Dept. for X-ray pad part numbers and additional information.



Astrus Elite Silver™ OR Table Pad Sets

OR Table Manufacturer / Model Name	Astrus Elite Silver - Order Info			Astrus Elite Silver Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
Berchtold B810, B710	H-30A-12-100	H-40A-12-100	H-50A-12-100	H-31A-12-100	H-41A-12-100	H-51A-12-100
Berchtold Operon D850, D750 (long foot)	H-30A-12-120	H-40A-12-120	H-50A-12-120	H-31A-12-120	H-41A-12-120	H-51A-12-120
Berchtold Operon D850, D750 (short foot)	H-30A-12-130	H-40A-12-130	H-50A-12-130	H-31A-12-130	H-41A-12-130	H-51A-12-130
Maquet 1130	H-30A-10-110	H-40A-10-110	H-50A-10-110	H-31A-10-110	H-41A-10-110	H-51A-10-110
Maquet 1130 Split Leg	H-30A-10-100	H-40A-10-100	H-50A-10-100	H-31A-10-100	H-41A-10-100	H-51A-10-100
Maquet 1130 Orthopedic	H-30A-10-130	H-40A-10-130	H-50A-10-130	H-31A-10-130	H-41A-10-130	H-51A-10-130
Maquet Alphamaxx (long torso)	H-30A-10-150	H-40A-10-150	H-50A-10-150	H-31A-10-150	H-41A-10-150	H-51A-10-150
Maquet Alphamaxx (short torso)	H-30A-10-160	H-40A-10-160	H-50A-10-160	H-31A-10-160	H-41A-10-160	H-51A-10-160
Maquet Alphastar 1132	H-30A-10-120	H-40A-10-120	H-50A-10-120	H-31A-10-120	H-41A-10-120	H-51A-10-120
Midmark 7100	H-30A-06-100	H-40A-06-100	H-50A-06-100	H-31A-06-100	H-41A-06-100	H-51A-06-100
Midmark Chic CST	H-30A-06-120	H-40A-06-120	H-50A-06-120	H-31A-06-120	H-41A-06-120	H-51A-06-120
Midmark Montreal	H-30A-06-130	H-40A-06-130	H-50A-06-130	H-31A-06-130	H-41A-06-130	H-51A-06-130
Getinge / Champagne 4900B, 4900E	H-30A-04-100	H-40A-04-100	H-50A-04-100	H-31A-04-100	H-41A-04-100	H-51A-04-100
Getinge / Champagne 5100B, 5100E	H-30A-04-130	H-40A-04-130	H-50A-04-130	H-31A-04-130	H-41A-04-130	H-51A-04-130
Siemens Universal	H-30A-13-100	H-40A-13-100	H-50A-13-100	H-31A-13-100	H-41A-13-100	H-51A-13-100
Skytron 5001, 6001, 6500, 6600	H-30A-05-130	H-40A-05-130	H-50A-05-130	H-31A-05-130	H-41A-05-130	H-51A-05-130
Skytron 3100	H-30A-05-100	H-40A-05-100	H-50A-05-100	H-31A-05-100	H-41A-05-100	H-51A-05-100
Skytron 3500	H-30A-05-120	H-40A-05-120	H-50A-05-120	H-31A-05-120	H-41A-05-120	H-51A-05-120
Skytron 3600	H-30A-05-190	H-40A-05-190	H-50A-05-190	H-31A-05-190	H-41A-05-190	H-51A-05-190
Skytron 6700, 6701	H-30A-05-200	H-40A-05-200	H-50A-05-200	H-31A-05-200	H-41A-05-200	H-51A-05-200
Steris-Amsco 1080, 2080	H-30A-03-100	H-40A-03-100	H-50A-03-100	H-31A-03-100	H-41A-03-100	H-51A-03-100
Steris-Amsco 3080, 3085	H-30A-03-110	H-40A-03-110	H-50A-03-110	H-31A-03-110	H-41A-03-110	H-51A-03-110
Steris Cmax / Steris 4085	H-30A-03-150	H-40A-03-150	H-50A-03-150	H-31A-03-150	H-41A-03-150	H-51A-03-150
Steris SurgiGraphic 6000	H-30A-03-170	H-40A-03-170	H-50A-03-170	H-31A-03-170	H-41A-03-170	H-51A-03-170

Individual pads for these sets are available. Contact Sales Dept. for X-ray pad part numbers and additional information.



Standard & Deluxe Surfaces

Deluxe Surfaces: Our conductive Deluxe Surface features supportive foam and a stretchable cover for added patient comfort. Recommended for use during shorter procedures and with patients who are determined to be at lower risk for developing pressure ulcers. Velcro® kit included.


Standard Surfaces: Our conductive Standard Surface is equivalent to the pads that are typically provided with new surgical tables. It combines supportive urethane foam with a 'Lectrolite® cover.

Standard & Deluxe OR Table Pad Sets


	Standard - Order Info		Deluxe - Order Info		Deluxe Gel - Order Info	
	2" Pad	3" Pad	2" Pad	3" Pad	2.25" Pad	3.25" Pad
OR Table Manufacturer / Model Name	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Berchtold B810, B710	H-201-12-100	H-301-12-100	H-202-12-100	H-302-12-100	H-212-12-100	H-312-12-100
Berchtold Operon D850, D750 (long foot)	H-201-12-120	H-301-12-120	H-202-12-120	H-302-12-120	H-212-12-120	H-312-12-120
Berchtold Operon D850, D750 (short foot)	H-201-12-130	H-301-12-130	H-202-12-130	H-302-12-130	H-212-12-130	H-312-12-130
Maquet 1130	H-201-10-110	H-301-10-110	H-202-10-110	H-302-10-110	H-212-10-110	H-312-10-110
Maquet 1130 Split Leg Table	H-201-10-100	H-301-10-100	H-202-10-100	H-302-10-100	H-212-10-100	H-312-10-100
Maquet 1130 Orthopedic Table	H-201-10-130	H-301-10-130	H-202-10-130	H-302-10-130	H-212-10-130	H-312-10-130
Maquet Alphamaxx (long torso)	H-201-10-150	H-301-10-150	H-202-10-150	H-302-10-150	H-212-10-150	H-312-10-150
Maquet Alphamaxx (short torso)	H-201-10-160	H-301-10-160	H-202-10-160	H-302-10-160	H-212-10-160	H-312-10-160
Maquet Alphastar 1132	H-201-10-120	H-301-10-120	H-202-10-120	H-302-10-120	H-212-10-120	H-312-10-120
Midmark 7100	H-201-06-100	H-301-06-100	H-202-06-100	H-302-06-100	H-212-06-100	H-312-06-100
Midmark Chic CST	H-201-06-120	H-301-06-120	H-202-06-120	H-302-06-120	H-212-06-120	H-312-06-120
Midmark Montreal	H-201-06-130	H-301-06-130	H-202-06-130	H-302-06-130	H-212-06-130	H-312-06-130
Getinge / Champagne 4900B, 4900E	H-201-04-100	H-301-04-100	H-202-04-100	H-302-04-100	H-212-04-100	H-312-04-100
Getinge / Champagne 5100B, 5100E	H-201-04-130	H-301-04-130	H-202-04-130	H-302-04-130	H-212-04-130	H-312-04-100
Siemens Universal	H-201-13-100	H-301-13-100	H-202-13-100	H-302-13-100	H-212-13-100	H-312-13-100
Skytron 5001, 6001, 6500, 6600	H-201-05-130	H-301-05-130	H-202-05-130	H-302-05-130	H-212-05-130	H-312-05-130
Skytron 3100	H-201-05-100	H-301-05-100	H-202-05-100	H-302-05-100	H-212-05-100	H-312-05-100
Skytron 3500	H-201-05-120	H-301-05-120	H-202-05-120	H-302-05-120	H-212-05-120	H-312-05-120
Skytron 3600	H-201-05-190	H-301-05-190	H-202-05-190	H-302-05-190	H-212-05-190	H-312-05-190
Skytron 6700, 6701	H-201-05-200	H-301-05-200	H-202-05-200	H-302-05-200	H-212-05-200	H-312-05-200
Steris-Amsco 1080, 2080	H-201-03-100	H-301-03-100	H-202-03-100	H-302-03-100	H-212-03-100	H-312-03-100
Steris-Amsco 3080, 3085	H-201-03-110	H-301-03-110	H-202-03-110	H-302-03-110	H-212-03-110	H-312-03-110
Steris Cmax / Steris 4085	H-201-03-150	H-301-03-150	H-202-03-150	H-302-03-150	H-212-03-150	H-312-03-150
Steris SurgiGraphic 6000	H-201-03-170	H-301-03-170	H-202-03-170	H-302-03-170	H-212-03-170	H-312-03-170

Note: Pad products, including gel pads, FlexiForm™, Allen® Hug-U-Vac® and foam positioners, returned within thirty (30) days of the invoice date are subject to a 20% restocking fee.


OR Table Pad Sizing Chart




Berchtold B810, B710




Maquet 1130




Berchtold D850, D750 (long foot)




Maquet 1130 Split Leg




Berchtold D850, D750 (short foot)




Maquet 1130 Orthopedic



Getinge / Champagne 4900B




Maquet Alphastar 1132



Maquet Alphamaxx (long torso)



Midmark 7100




Maquet Alphamaxx (short torso)




Midmark Chic CST


OR Table Pad Sizing Chart




Midmark Montreal




Skytron 5001 / 6001 / 6500 / 6600




Siemens Universal




Skytron 6700




Skytron 3500




Steris/Amsco 1080 / 2080




Steris/Amsco 3080 / 3085




Skytron 3600




Steris Cmax / Steris 4085




Choose Your Armboard Style



Fits armboard models:
• Allen® Japanese



Fits most armboard models:
• Getinge® • MDT
• Maquet® • Shampaine®



Fits most armboard models:
• Allen • Schaerer Mayfield®
• Amsco® • Siemens®
• Berchtold® • Steris®



Fits armboard models:
• Skytron®

Call Sales Dept. for information on additional armboard model pads.
All armboard pads sold individually.



Alto™ Armboards Surfaces

Effective Pressure Management

Engineered, multi-layered armboard surfaces are designed to reduce pressure and provide support for the limb. This support system decreases the high pressure that bony prominences usually experience. The anti-shear cover helps protect patient's skin.

Allen® armboard surfaces are available for armboards in both regular and contoured styles. Allen contoured armboard surface rounds up at the sides an extra ½" for increased stability.

All armboard pads sold individually.

Standard Surfaces



	Alto Order Info		Alto Gel Order Info	
	4" Pad	5" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-006-03-950	-----	H-016-03-950	-----
Berchtold®	H-006-12-950	-----	H-016-12-950	-----
Maquet®	H-006-10-950	-----	H-016-10-950	-----
Midmark®	H-006-06-950	-----	H-016-06-950	-----
Getinge™ / Champagne®	H-006-04-950	-----	H-016-04-950	-----
Skytron®	H-006-05-950	-----	H-016-05-950	-----
Steris-Amsco X-ray	-----	H-006-03-951	-----	H-016-03-951
Skytron X-ray	-----	H-006-05-951	-----	H-016-05-951

Contoured Surfaces



	Alto Order Info		Alto Gel Order Info	
	4" Pad	5" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-006-03-960	-----	H-016-03-960	-----
Berchtold	H-006-12-960	-----	H-016-12-960	-----
Maquet	H-006-10-960	-----	H-016-10-960	-----
Midmark	H-006-06-960	-----	H-016-06-960	-----
Getinge / Champagne	H-006-04-960	-----	H-016-04-960	-----
Skytron	H-006-05-960	-----	H-016-05-960	-----
Steris-Amsco X-ray	-----	H-006-03-961	-----	H-016-03-961
Skytron X-ray	-----	H-006-05-961	-----	H-016-05-961

Note: "-----" means pad is not available.



Astrus Elite™ & Astrus Elite Silver™ Armboards Surfaces

Effective Pressure Management - Astrus Elite and Astrus Elite Silver Surfaces have a 4-layer foam composite construction that helps protect patients four ways. Allen® armboard surfaces are available for armboards in both regular and contoured styles. Allen contoured armboard surface rounds up at the sides an extra ½" for increased stability.

All armboard pads sold individually.

Standard Surfaces



	Astrus Elite - Order Info			Astrus Elite Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-303-03-950	H-403-03-950	H-503-03-950	H-313-03-950	H-413-03-950	H-513-03-950
Berchtold	H-303-12-950	H-403-12-950	H-503-12-950	H-313-12-950	H-413-12-950	H-513-12-950
Maquet	H-303-10-950	H-403-10-950	H-503-10-950	H-313-10-950	H-413-10-950	H-513-10-950
Midmark	H-303-06-950	H-403-06-950	H-503-06-950	H-313-06-950	H-413-06-950	H-513-06-950
Getinge / Shampaine	H-303-04-950	H-403-04-950	H-503-04-950	H-313-04-950	H-413-04-950	H-513-04-950
Skytron	H-303-05-950	H-403-05-950	H-503-05-950	H-313-05-950	H-413-05-950	H-513-05-950

Contoured Surfaces



	Astrus Elite - Order Info			Astrus Elite Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-303-03-960	H-403-03-960	H-503-03-960	H-313-03-960	H-413-03-960	H-513-03-960
Berchtold	H-303-12-960	H-403-12-960	H-503-12-960	H-313-12-960	H-413-12-960	H-513-12-960
Maquet	H-303-10-960	H-403-10-960	H-503-10-960	H-313-10-960	H-413-10-960	H-513-10-960
Midmark	H-303-06-960	H-403-06-960	H-503-06-960	H-313-06-960	H-413-06-960	H-513-06-960
Getinge / Shampaine	H-303-04-960	H-403-04-960	H-503-04-960	H-313-04-960	H-413-04-960	H-513-04-960
Skytron	H-303-05-960	H-403-05-960	H-503-05-960	H-313-05-960	H-413-05-960	H-513-05-960

Standard Surfaces



	Astrus Elite Silver - Order Info			Astrus Elite Silver Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-30A-03-950	H-40A-03-950	H-50A-03-950	H-31A-03-950	H-41A-03-950	H-51A-03-950
Berchtold	H-30A-12-950	H-40A-12-950	H-50A-12-950	H-31A-12-950	H-41A-12-950	H-51A-12-950
Maquet	H-30A-10-950	H-40A-10-950	H-50A-10-950	H-31A-10-950	H-41A-10-950	H-51A-10-950
Midmark	H-30A-06-950	H-40A-06-950	H-50A-06-950	H-31A-06-950	H-41A-06-950	H-51A-06-950
Getinge / Shampaine	H-30A-04-950	H-40A-04-950	H-50A-04-950	H-31A-04-950	H-41A-04-950	H-51A-04-950
Skytron	H-30A-05-950	H-40A-05-950	H-50A-05-950	H-31A-05-950	H-41A-05-950	H-51A-05-950

Contoured Surfaces



	Astrus Elite Silver - Order Info			Astrus Elite Silver Gel - Order Info		
	3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-30A-03-960	H-40A-03-960	H-50A-03-960	H-31A-03-960	H-41A-03-960	H-51A-03-960
Berchtold	H-30A-12-960	H-40A-12-960	H-50A-12-960	H-31A-12-960	H-41A-12-960	H-51A-12-960
Maquet	H-30A-10-960	H-40A-10-960	H-50A-10-960	H-31A-10-960	H-41A-10-960	H-51A-10-960
Midmark	H-30A-06-960	H-40A-06-960	H-50A-06-960	H-31A-06-960	H-41A-06-960	H-51A-06-960
Getinge / Shampaine	H-30A-04-960	H-40A-04-960	H-50A-04-960	H-31A-04-960	H-41A-04-960	H-51A-04-960
Skytron	H-30A-05-960	H-40A-05-960	H-50A-05-960	H-31A-05-960	H-41A-05-960	H-51A-05-960



Standard & Deluxe Armboards Surfaces

Deluxe Surfaces: Recommended for use during shorter procedures and with patients who are determined to be at lower risk for developing pressure ulcers. Velcro® kit included.

Standard Surfaces: Our conductive Standard Surface is equivalent to the pads that are typically provided with new surgical tables.

All armboard pads sold individually.



Standard Surfaces

	Standard - Order Info			Deluxe - Order Info			Deluxe Gel - Order Info		
	2" Pad	3" Pad	4" Pad	2" Pad	3" Pad	4" Pad	2" Pad	3" Pad	4" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen®, Steris-Amsco	H-201-03-950	H-301-03-950	H-401-03-950	H-202-03-950	H-302-03-950	H-402-03-950	H-212-03-950	H-312-03-950	H-412-03-950
Berchtold	H-201-12-950	H-301-12-950	H-401-12-950	H-202-12-950	H-302-12-950	H-402-12-950	H-212-12-950	H-312-12-950	H-412-12-950
Maquet	H-201-10-950	H-301-10-950	H-401-10-950	H-202-10-950	H-302-10-950	H-402-10-950	H-212-10-950	H-312-10-950	H-412-10-950
Midmark	H-201-06-950	H-301-06-950	H-401-06-950	H-202-06-950	H-302-06-950	H-402-06-950	H-212-06-950	H-312-06-950	H-412-06-950
Getinge / Champagne	H-201-04-950	H-301-04-950	H-401-04-950	H-202-04-950	H-302-04-950	H-402-04-950	H-212-04-950	H-312-04-950	H-412-04-950
Skytron	H-201-05-950	H-301-05-950	H-401-05-950	H-202-05-950	H-302-05-950	H-402-05-950	H-212-05-950	H-312-05-950	H-412-05-950



Contoured Surfaces

	Standard - Order Info			Deluxe - Order Info			Deluxe Gel - Order Info		
	2" Pad	3" Pad	4" Pad	2" Pad	3" Pad	4" Pad	2" Pad	3" Pad	4" Pad
	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	Product No.
Allen, Steris-Amsco	H-201-03-960	H-301-03-960	H-401-03-960	H-202-03-960	H-302-03-960	H-402-03-960	H-212-03-960	H-312-03-960	H-412-03-960
Berchtold	H-201-12-960	H-301-12-960	H-401-12-960	H-202-12-960	H-302-12-960	H-402-12-960	H-212-12-960	H-312-12-960	H-412-12-960
Maquet	H-201-10-960	H-301-10-960	H-401-10-960	H-202-10-960	H-302-10-960	H-402-10-960	H-212-10-960	H-312-10-960	H-412-10-960
Midmark	H-201-06-960	H-301-06-960	H-401-06-960	H-202-06-960	H-302-06-960	H-402-06-960	H-212-06-960	H-312-06-960	H-412-06-960
Getinge / Champagne	H-201-04-960	H-301-04-960	H-401-04-960	H-202-04-960	H-302-04-960	H-402-04-960	H-212-04-960	H-312-04-960	H-412-04-960
Skytron	H-201-05-960	H-301-05-960	H-401-05-960	H-202-05-960	H-302-05-960	H-402-05-960	H-212-05-960	H-312-05-960	H-412-05-960



Astrus Elite™ & Astrus Elite Silver™ Stretcher Surfaces

Effective Pressure Management 3½" (9 cm) Surface

Our Astrus Elite Surface has a 4-layer foam composite construction that helps protect patients four ways. A top layer of memory foam conforms to the patient while the underlying foam layers support the patient's weight and prevent bottoming out. Both the Astrus Elite and Astrus Elite Silver contain the same 4-layer advanced pad technology.

Astrus Elite Stretcher Surfaces

Width	Length: 74" (188 cm)	Length: 76" (193 cm)	Length: 78" (198 cm)
	3.5" Pad	3.5" Pad	3.5" Pad
	Product Numbers	Product Numbers	Product Numbers
24"	H-323-00-424	H-323-00-624	H-323-00-824
26"	H-323-00-426	H-323-00-626	H-323-00-826
28"	H-323-00-428	H-323-00-628	H-323-00-828
30"	H-323-00-430	H-323-00-630	H-323-00-830

Astrus Elite Silver Stretcher Surfaces

Width	Length: 74" (188 cm)	Length: 76" (193 cm)	Length: 78" (198 cm)
	3.5" Pad	3.5" Pad	3.5" Pad
	Product Numbers	Product Numbers	Product Numbers
24"	H-32A-00-424	H-32A-00-624	H-32A-00-824
26"	H-32A-00-426	H-32A-00-626	H-32A-00-826
28"	H-32A-00-428	H-32A-00-628	H-32A-00-828
30"	H-32A-00-430	H-32A-00-630	H-32A-00-830

Individual pads for these sets are available. Contact Sales Dept. for X-ray pad part numbers and additional information.



Standard & Deluxe Stretcher Surfaces

Cost Effective 2" (5 cm) Surface 3½" (9 cm) Surface

The Standard Pad features supportive foam and a 'Lectrolite®' cover. This surface is equivalent to standard pads that are usually provided with stretchers.

The Deluxe Pad features supportive foam and a stretchable cover for added comfort.

Deluxe Stretcher Surfaces

Width	Length: 74" (188 cm)		Length: 76" (193 cm)		Length: 78" (198 cm)	
	2" Pad	3.5" Pad	2" Pad	3.5" Pad	2" Pad	3.5" Pad
	Product Numbers		Product Numbers		Product Numbers	
24"	H-222-00-424	H-322-00-424	H-222-00-624	H-322-00-624	H-222-00-824	H-322-00-824
26"	H-222-00-426	H-322-00-426	H-222-00-626	H-322-00-626	H-222-00-826	H-322-00-826
28"	H-222-00-428	H-322-00-428	H-222-00-628	H-322-00-628	H-222-00-828	H-322-00-828
30"	H-222-00-430	H-322-00-430	H-222-00-630	H-322-00-630	H-222-00-830	H-322-00-830

Standard Stretcher Surfaces

Width	Length: 74" (188 cm)		Length: 76" (193 cm)		Length: 78" (198 cm)	
	2"	3.5"	2"	3.5"	2"	3.5"
	Product Numbers		Product Numbers		Product Numbers	
24"	H-221-00-424	H-321-00-424	H-221-00-624	H-321-00-624	H-221-00-824	H-321-00-824
26"	H-221-00-426	H-321-00-426	H-221-00-626	H-321-00-626	H-221-00-826	H-321-00-826
28"	H-221-00-428	H-321-00-428	H-221-00-628	H-321-00-628	H-221-00-828	H-321-00-828
30"	H-221-00-430	H-321-00-430	H-221-00-630	H-321-00-630	H-221-00-830	H-321-00-830

Allen® Gel Positioners

Features

- Constructed with an advanced silicone formula designed to help with pressure redistribution
- Flexible membrane cover
- Will not absorb body fluids or odors
- Radiolucent

Armboard & Toboggan Pads



Small Armboard Gel Pad
15½" x 4½" x ¾"
(395 x 115 x 10 mm)
#A-10105



Large Armboard Gel Pad
23¾" x 5" x ¾"
(600 x 125 x 10 mm)
#A-10032



Allen® Armboard Gel Pad
23¾" x 8" x ¾"
(600 x 200 x 10 mm)
#A-10035



Allen Premium Armboard Gel Pad - Short
13" x 5½" x 1.2"
(330 x 140 x 30 mm)
#A-10036



Allen Premium Armboard Gel Pad - Long
18½" x 7" x 1½"
(460 x 180 x 35 mm)
#A-10037



Toboggan Gel Pad
15½" x 10" x ¾"
(395 x 245 x 10 mm)
#A-10029

Head Positioners



Adult Head Gel Donut
8" Outer Diameter x 2"
(200 x 45 mm)
#A-10031



Pediatric Head Gel Donut
3½" Outer Diameter x ¾"
(90 x 20 mm)
#A-10103



Adult Gel Horseshoe
8" Outer Diameter x 3"
(200 x 75 mm)
#A-10104



Allen® Prone Head Positioner - Large
11" x 9½" x 5½"
(280 x 240 x 140 mm)
#A-10038



Allen Prone Head Positioner - Small
11" x 9" x 5½"
(273 x 220 x 135 mm)
#A-10039



Allen Silicone Head Positioner
8" x 10"
(203 x 254 mm)
#A-10045

Knee Pads *Continued On Page 140*

Small Knee Crutch Gel Pad
12⅝" x 8⅞" x ¾"
(320 x 220 x 10 mm)
#A-10070



Medium Knee Crutch Gel Pad
15½" x 10" x ¾"
(395 x 245 x 10 mm)
#A-10029



Hip/Large Knee Crutch Gel Pad
20½" x 20½" x ¾"
(520 x 520 x 10 mm)
#A-10033

Knee Pads *Continued*



Allen® Prone Knee Pad - Long
20½" x 6" x 1½"
(510 x 150 x 35 mm)
#A-10040



Allen Prone Knee Pad - Short
16" x 6" x 1½"
(400 x 150 x 35 mm)
#A-10041

Nerve & Heel Protectors



Small Ulnar Nerve Gel Elbow Pad With Velcro® Strap
15" x 4½" x 5/16"
(380 x 115 x 8 mm)
#A-10102



Large Ulnar Nerve Gel Elbow Pad With Velcro Straps
18" x 8¼" x 5/16"
(460 x 210 x 8 mm)
#A-10066



Foot-Heel Gel Pad With Velcro Straps
3/8" x 3¾" x 4½"
(10 x 95 x 114 mm)
#A-10067



Heel Protector Gel Pad With Velcro Straps
9" x 7⅞" x 7/8"
(230 x 200 x 20 mm)
#A-10068



Polymer Gel Boot Pads
19⅞" x 5¾" x 3/8"
(505 x 145 x 10 mm)
#A-10026

Positioners



“Sandbag” Gel Positioner
Small: 12" x 4¾" x 2¾" (300 x 120 x 60 mm),
Large: 14" x 5" x 2¾" (360 x 125 x 70 mm)
#A-10073 (Small)
#A-10074 (Large)



Chest Roll Gel Positioner
12" x 4" x 3" (305 x 100 x 75 mm), 16" x 4" x 3"
(405 x 100 x 75 mm), 20" x 4" x 3" (510 x 100 x 75 mm)
#A-10076
#A-10078
#A-10080

Positioners *Continued*


Chest Roll Gel Positioner With Foam Core
 16" x 6" x 6" (405 x 150 x 150 mm),
 20" x 6" x 6" (510 x 150 x 150 mm)
 #A-10097
 #A-10098



Allen® Prone Gel Positioner
 44½" x 21" x 5½"
 (1115 x 520 x 140 mm)
 #A-10044



Allen Lateral Gel Positioner
 29" x 20" x 6½"
 (720 x 500 x 160 mm)
 #A-10043



Allen Lateral Leg Positioner
 26" x 16" x 9"
 (650 x 400 x 230 mm)
 #A-10042

Overlay Table Pads




Universal 3-Piece OR Table
 Gel Pad Set
 #A-10034



Head Section Gel Pad
 9" x 20" x ¾"
 (225 x 500 x 10 mm)
 #A-10034H



Torso Section Gel Pad With Cutout
 39" x 20" x ¾"
 (990 x 500 x 10 mm)
 #A-10034T



Foot Section Gel Pad
 19" x 20" x ¾"
 (500 x 500 x 10 mm)
 #A-10034F




¾ Length: 45¼" x 20½" x ¾" (1150 x 520 x 10 mm),
 Full Length: 71" x 20½" x ¾" (1800 x 520 x 10 mm)
 #A-10060 (¾ Length) Table Gel Pad
 #A-10061 (Full Length) Table Gel Pad

Lateral Positioning



- A** #A-10043 Allen® Lateral Gel Positioner
- B** #A-10066 Large Ulnar Nerve Gel Elbow Pad w/ Velcro® straps



- C** #A-10037 Allen Premium Armboard Gel Pad - Long
- D** #A-10029 Toboggan Gel Pad
- E** #A-10042 Allen Lateral Leg Positioner

Prone Positioning



- A** #A-10040 Allen Prone Knee Pad - Long
- B** #A-10097 Chest Roll Gel Positioner w/ Foam Core
- C** #A-10044 Allen Prone Gel Positioner



- D** #A-10038 Allen Prone Head Positioner - Large
- E** #A-10037 Allen Premium Armboard Gel Pad - Long