

# Surgical Stretcher - P8010

Efficient surgical patient transport and positioning



## Features and benefits

**Heavy-gauge, steel-reinforced frame** ensures durability, even with frequent usage

**Fully integrated siderails** mean safer patient transport

**Articulating head section** for maximum flexibility

**Surgical rail** for easy attachment of gas, extenders and hand tables

**Steering Plus™** directional system incorporates fully-locking five castor braking, providing improved steering, stopping, cornering and control

**Integrated handbrake system** offers power-assisted braking

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**


## Surgical stretcher - P8010



Optional: push handles



Articulating head section



Armboard extension

### Technical Specifications (standard configuration)

Safe Working Load	317.5 kg
Weight (without surface/accessories)	129.3 kg
Height adjustment	54.6 to 84.5 cm
Trendelenburg/Reverse Trendelenburg	- 18 to + 18 °
Protection class	IPX4
Deck width	66 cm
<b>Width</b>	
siderails stored	74.3 cm
siderails up	81.3 cm
Deck length	198.1 cm
Length	233.7 cm
Surface height	10.2 cm
<b>Siderails</b>	
height	36.8 cm
length	119.4 cm

Clearance	2.9 to 8.9 cm
Castors	4 x 20.3 cm single band

### Additional options

Push handles at foot end only  
Hydraulic control pedals on both sides at foot end  
4 brake and steer pedals

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This product is not available in countries requiring a CE mark. This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN145304-02, 22 March 2016

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**